


29. Radial Division of an Object

Like grid division, radial division of objects can be useful for subdividing structures into various substructures.

- Select MRI_3D_Head and open Segment.
- Select Manual **1** and choose Auto Trace. **2**
- Set the primary display to Axial **3** and use Slice **4** to move to axial slice 130.
- Set a seed pixel in the white matter **5** and adjust the minimum and maximum threshold values **6** to define the brain. A min value of 45 and max value of 125 works well for this data set.
- Click Apply **7** to define the brain.



Download the MRI_3D_Head data set to follow along
<http://analyzedirect.com/data/>



- 10

