

## 6. Making Angle Measurements

The Angles tool allows users to make distance measurements between points and angle measurements between defined lines.



Download the MRI\_3D\_Head data set to follow along <http://analyzedirect.com/data/>

- Select MRI\_3D\_Head and open Measure. **1**
- Select File > Load Object Maps and load MRI\_3D\_Head.obj.
- Select the 1D Sample Type **1** and choose Angles. **2**
- Check the Points, Distances, and Angles to be reported. **3**
- Define the angles on any of the 2D slices or the 3D rendering. **4**
- The selected measurements will be reported in the Stats review area. **5**

